



— Pharma Process & Packaging

www.saintyco.com

EcoTab-15Plus

Rotary Tablet Press



EcoTab-15Plus

Rotary tablet press



EASY
TO USE

MAX.
TABLET
50MM

UP TO
9000
TABLETS

EcoTab-15Plus rotary tablet press is used to press all sorts of granular materials, like geometry shape of profiled sheet, double lettering piece, traditional Chinese medicine tablet, western medicine tablet, calcium tablet, buccal tablet, spirulina tablet, milk tablet, candy tablet, effervescent tablets, chewable tablets, seasoning block, chicken pieces, camphor ball (naphthalene), aromatic piece, disinfection tablet, fertilizer tablet, pesticide tablet, button cell, electronics, metallurgy powder ceramic grain, etc.

Features



Features

- This machine is designed for pressing round tablets from various kinds of granular materials. It is primary machine to fit batch production of tablets.
- During operation, the speed of the turret rotation, the depth of the filling material and the thickness of the tablets can be adjusted.
- The mechanical buffer unit will prevent punches and apparatus from overload damage.
- A powder suction unit is attached to take in the residual powder can be eliminated and the collected powder material can be reused.

Technical Data



Technical Data	EcoTab-15Plus
Quantity of Punch & Die (set)	15
Max. Pressure (kn)	120
Max. Dia. of Tablet (mm)	50
Max. Thickness of Tablet (mm)	18
Max. Depth of Filling (mm)	35
Rotate Speed (r/min)	0-20
Production Capacity (pcs/h)	38400
Motor Power (kw)	7.5
Overall Size (mm)	1250×980×1650
Weight (kg)	1850

SaintyCo also offer:

- Fluid-Bed Processors
- High Shear Mixer Granulator
- Roller Compactor
- Blister Packing Solutions

For additional information, please drop email to info@saintyco.com.

