

Drum Lifting Column-Fixed only in the floor

NTG Series



Drum Lifting Column-Fixed only in the floor

NTG series Drum Lifting Column is mainly used in the solid dosage production procedure for transferring and discharging materials, consequently it is operated with corollary equipments such as granulator, blender, tablets press, capsule filling machine, etc. Pharma Lifters, stationary are widely used in pharmaceutical, chemical, foodstuff industries, etc.




Features:

- Max load: 300 kg.
- It composed of chassis, column and drum inversion system.
- It has a reasonable structure with stable performances and easy operation.
- There has no dead angle or protruded screws.
- It designed with a device avoiding misoperation of the butterfly valve.
- Easy to wash.
- cGMP compliant.



Dimensions:



Technical Data												
Model	Main structural dimensions									Lifting Capacity (with drum)	Total power	Machine weight
	L	L1	L2	W	H	H1	H2	H3	DN	kg	kw	kg
NTG50	1800	420	110 0	60 0	320 0	220 0	255 0	298 5	15 0	50	1.5	270
NTG100	1920	500	110 0	70 0	340 0	220 0	255 0	315 0	20 0	100	2	340
NTG150	1920	500	110 0	70 0	351 0	220 0	255 0	326 0	20 0	150	2.2	340
NTG200	1920	500	110 0	70 0	351 0	220 0	255 0	332 0	20 0	200	2.5	340
NTG300	1920	500	110 0	70 0	315 0	220 0	255 0	350 0	20 0	300	2.5	340

SaintyCo also offer:

Encapsulation Solutions

Tablet Compression Solutions

Coating Systems

Cartoning Systems

Blister Packing Solutions

